



METAL-FREE EARLOOP MASK

Item No. FM81001



MASK TYPE	Procedure Earloop
MASK COLOR	Tan
MASK MATERIAL	Cellulose
EARLOOP MATERIAL	Elastic
NOSEPIECE	Elastic
LATEX CONTENT	Not made with natural rubber latex
METAL CONTENT	Does not contain metal
DELTA P	2.2
PFE	+98%
BFE	+98%
ITEMS PER BOX	40 masks
BOXES PER CASE	10 boxes
ITEMS PER CASE	400 masks
CASE DIMENSIONS	19½" x 8½" x 8"
CASE WEIGHT	6 lbs
COUNTRY OF ORIGIN	USA